Experiment 9: Direct Titration of Lead with ErioT and EDTA

- Synopsis Lead is held in solution by weakly chelating tartrate so that it may react with Erio T to form a bluish violet color. The lead is titrated with standard EDTA. The end point is observed by a loss of the bluish violet color as the last of the lead-Erio-T complex is consumed.
- **<u>READINGS</u>** Read pages 279-285 in Critical Reviews

Reagents

- 0.01 M lead solution
- 0.01 M EDTA standard solution: Dilute from 0.100 M EDTA
 - 37.22 g of disodium EDTA in 1000 mL deionest water.

Or

Dry the acid for two hours at 130-150C. Cool. Weigh 29.210 g of acid EDTA, add to 600 mL water, add pellet by pellet NaOH, until the EDTA comes into solution. Dilute to 1L.

Erio T indicator powder:

Grind 100 mg of indicator with 10 g of NaCl to a very fine powder and store.

Tartaric acid

pH 10 buffer

: Dissolve 70 g of NH_4Cl in 570 ml of ammonia (s. G. 0.90) and dilute to distilled water to 1 L.

Procedure

- 1. Place 10-30 ml (exactly measured) 0.01 M lead soluiton in 250 ml flask
- 2. Add a spatula end of tartaric acid.
- 3. Add 5 ml of buffer pH 10 and dilute to about 50-100 ml. If a turbidity occurs (Pb(OH)₂) add more tartaric acid.
- 4. Add Erio T (too much will change the color intensity, so start small).
- 5. Titrate until the colour changes from violet just to clear blue.
- 6. Repeat twice to be able to report the rsd of the method.

1 ml 0.01 M EDTA = 2.0719 mg Pb

<u>REPORT</u> In addition to material, methods and results, include:

- 1. What is the rsd of this method?
- 2. What will determine the minimum amount of lead that can be measured in this method?
- 3. What constitutes a blank in this procedure? What are the sources of error embodied in the

standard deviation of the blank?

- 4. How does sample matrix affect your results?
- 6. What was the estimated time for turn around in samples?
- 7. Are there any problems with disposal of hazardous materials?
- 8. How easy would it be to instruct a technician on this method?
- 9. How easy would it be to construct a paper trail for this method?
- 10. What metals will interfere the most with Pb in this method?
- 11. Why doesn't the lead precipitate. Use the actual values of complex formation constants to justify your answer.
- 12. What is the color forming species and why is it affected by lead?
- 13. Why does the pH have to be brought to 12?